

Down-conversion crystals (double crystal)	Newlight Photonics	Paired BBO (2pcs) for photon entanglement, size 5x5x0.5 mm(each), cut for Type I PDC pumped by 405nm with the external half opening angle of 3 deg, theta=29.3deg., AR coated at 810/405 nm, mounted back-to-back with one crystal rotated by 90 degree about axis normal to incident face in an 1" anodized Al holder	PABBO5050-405(I)-HA3	1	1079	1079	
0.5mm BBO crystal rochon polarizer	Newlight Photonics	x-cut 10x10x0.5 mm quartz substrate polished on 2 sides NOTE: for better results, replace this with thin BBO (optional equipment, below)	NCBBO5050-405(I)-HARPB0110	1	559	559	
polarizer mounts	Thorlabs	KM100		2	40	80	
and phase plate mount	Thorlabs	KM100		1	40	40	
crystal mount	Thorlabs	KM100		1	40	40	
base	Thorlabs	base	BA1	4	6	24	
2" post	Thorlabs	2" post	TR2	4	5	20	
2" post holder	Thorlabs	2" post holder	PH2	4	8	32	
Quantum State Measurement							
Quarter-wave plates	Newlightphotronics.com		WPA03-Q-810	2	309	618	
rotation stage for quarter-wave plates	Thorlabs		RSP1D	2	95	190	
Tools							
spanner wrench	Thorlabs	for Aspheric optics	SPW301	1	15	15	
spanner wrench	Thorlabs	For SM1 lens tubes	SPW602	1	27	27	
spanner wrench	Thorlabs	For fiber coupling lens focus adjustment	SPW909	1	29	29	
fiber microscope	thorlabs	cleaning fibers	fs201	1	224	224	
fiber cleaner	thorlabs	Cleans fiber tips	fcc-7020	1	23	23	
SUB-TOTAL							
						26833	This cost includes the Digikey detectors, not the ALPhA detectors
Optional Equipment							
CCD Camera	maxmax.com	usb camera, IR sensitive	XNiteUSB2S-M	2	85	170	very useful for doing alignment. It lets you see the 800 nm alignment laser on your computer monitor
TAC	Ortec	Time-to-Amplitude Converter for observing timing of photon pairs (NOTE: you'll also need a NIM crate to plug this into--most labs have one lying around somewhere)					This is pretty old-school, you'll also need an MCA
TDC	s-fifteen.com	time-to digital converter	TDC1	1	1500	1500	Only \$1,500, but also only 2ns resolution. I've never used this , but it should work if all you want to do is show that photon pairs come at the same time. You could also use this as your coincidence counter.
Green LED safe-lights							
TOTAL						28503	