

# Surface Plasmon Resonance: Exciting Waves and Metal Films

Nicole Moehring, Will Hollinger, Andrew Klavekoske, Todd Zimmerman\*

\*University of Wisconsin – Stout, Department of Chemistry and Physics

PHYS 439, Advanced Physics Lab

## Introduction:

Surface Plasmon Resonance is a popular technique that takes advantage of wave behavior at an interface to detect and/or identify materials and molecules. The purpose of this project was to answer the question “how can SPR be used to characterize thin metal films?” There were two main components to the project allowing the team to get closer to answering this question:

1. Measure gold film thickness using atomic force microscopy
2. Develop an experiment and setup that can be used to quantify plasmon resonance in response to three variables: film thickness, type of metal, and source light intensity

## Surface Plasmon Resonance (SPR):

SPR is often used to sense the presence of biomarkers. The phenomena is caused by the interaction of two wave types which form along the boundary of two dielectric materials:

- Evanescent wave – forms when light is totally internally reflected and decays exponentially as distance from the interface increases
- Surface plasmons – forms when free electrons in the metal film are excited and oscillating




Figure 1: Illustration explaining SPR technique<sub>1</sub>

## Measuring Film Thickness with AFM:

Atomic Force Microscopy (AFM): method for gathering accurate topography data on samples with surface feature sizes on the nanoscale




Figure 2: Illustration explaining atomic force microscopy and how it is used to measure film thickness by scanning across a thin film “step”


## Procedure:

- Gold film is deposited on a glass microscope slide using sputter coater ‘auto-cycle’ (30 second cycle)
- A scratch in the gold film is made, exposing the glass substrate and creating a step representative of the film thickness
- AFM is used to measure the step height and estimate the coating thickness each time the sputter coater ‘auto-cycle’ is run

## Results:

Cycle 10x = 97.62 nm  
Cycle 6x = 49.89 nm  
Cycle 2x = 25.41 nm

Figure 3: Boxplot of measured film thickness



## Quantifying Resonance:





Figure 3: Schematic of setup used to quantify plasmon resonance where L1/L2 are lenses, I1 is an iris, and PD is the photodetector

## Procedure:

- Deposit thin gold film on glass slide and place on edge of prism
- A red LED (634 nm) is used as the incident light source is shined through the prism where it interacts with the electrons in the film
- As the incident angle changes, a photodetector collects the reflected light and reports it as voltage
- When the evanescent waves and surface plasmons have equal energies, the incident light dissipates resulting in a decrease in intensity
- A relationship between incident angle and reflectivity and how the curve is affected by film thickness, film type, and incident light intensity can be graphed:

## Results:



## Conclusions:

- The sputter coater lays ~10 nm film during each ‘auto-cycle’
- Film material type changes the curve shape and incident angle where resonance occurs → helpful for characterizing film
- The thickness of the film changes the intensity of the resonance
- An iris opening smaller than 5.00 mm will affect reflectivity

## Acknowledgements:

The research team would like to thank Dr. Todd Zimmerman for valuable discussions and teaching throughout the project, Emily Frick for technical assistance on imaging equipment and the UW-Stout Department of Chemistry and Physics for funding.

## References:

1. Presentation, Kolomenski, A.; Surface Plasmons, Part 2, September 2012