




Turbulence: Flow Analysis By Imaging Particles

Mike Schatz, Bala Suri, Jeff Tithof (Georgia Tech)
Xiang Wu (Shanghai Jiao Tong University)

2012 Hands-On School:
17-29 June 2012
Shanghai Jiao Tong University
Shanghai, People's Republic of China




Today's Agenda

- Motivation
- Setting Up the Experiment
- Collecting Data
- Processing and Analyzing Data in Matlab

Motivation

- Statistical Treatment of Turbulence
- “Exact Coherent Structures” Guide Turbulence


Phase Space
Visualization



- <http://chaosbook.org/tutorials/>


Motivation

- Previous work
 - Hof et al. (2004)



- Experimental evidence is scarce
- Learn computational modeling in Dr. Storey's session (Mathematical Modeling Session B)
- Unique opportunity to compare experiment and numerics
- Your work on 2D turbulence could provide new evidence!



Setting Up the Experiment



Setting Up the Experiment

- Using a **Constant Current Source** is Ideal
 - Constant driving parameter
 - Prevent LEDs from flashing


50mA Constant Current Battery Charger



- More info in “ConstantCurrentSupply.pdf”

Setting Up the Experiment

- To visualize the flow
 - Sprinkle some tracer particles
 - Light up the flow with LEDs
 - Place a box over setup (to make it dark inside)
- Open “Fire-i” on the laptop
 - Start a “Preview”
 - Should look like this image:
- We will use a technique called “PIV” to extract velocities from these images



Collecting Data

- Adjust the variable resistor so that you see laminar flow ($\sim 2\text{mA}$ for Kolmogorov, $\sim 10\text{mA}$ for Checkerboard)
- Make a new folder on the desktop
- In Fire-i:
 - Make sure the frame rate is set to 7.5
 - Go to Settings->Frame Capture Properties, then set the folder and number of frames to capture
 - To start taking data, check “Frame Capture” and click “Start”

Collecting Data

- To see interesting dynamics, collect data at the following current values:

Kolmogorov	Checkerboard
2mA	10mA
7mA	25mA
15mA	45mA

Processing and Analyzing Data in Matlab

- Open Matlab, change to the folder containing your images, and execute these scripts:

Matlab Scripts (Listed in Execution Order)	Purpose
<code>image_process('startnum', 'stopnum', h_cm, fps);</code>	Calls OSIV software to perform PIV.
<code>taverage;</code>	Time averages data to reduce noise.
<code>vfield_cgs;</code>	Converts velocities and lengths to cgs units.
<code>vortex_cgs;</code>	Calculates vorticity contour in cgs units.
<code>remove_vortex_noise;</code>	Removes noisy parts of vorticity to make data look smoother.
<code>make_video(jump, video_fps);</code>	Generates an AVI video of the flow.
<code>recurrence;</code>	Generates a recurrence plot of the velocity field time-series.

- Examples:
 `image_process('002', '900', 4.8, 7.5);`
 `make_video(8, 15);`

Additional Analysis

- Calculating Reynolds numbers

$$\text{Re} = \frac{L_0 \cdot V_{rms}}{\nu} \quad \text{where} \quad \nu \approx 0.01 \text{ cm}^2 / \text{s}$$

- Time averaging steady state data

- Calculating Rayleigh friction decay constant

